

PHILIPS

Diagnostic ECG

Supplies catalog

Supplies catalog 2020





PHILIPS DigiTrak XT

PHILIPS
DigiTrak XT



Table of contents

Patient cables	2
Patient data cables	6
Replacement lead sets	7
Extended leads	8
Electrode connectors	9
Disposable ECG electrodes	10
Limb plate electrodes	11
Recording paper	11
Pouches and cases	13
Cardiograph batteries	14
Imaging supplies	15
TEE supplies	16
List of host products	17

Patient cables

Ordering code	Description	Packaging	Use with	CE	FDA
989803151631	Complete lead set, AAMI (2 arm leads, 2 leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 2 Trident yoke labels (V1-V3, V4-V6), 4 colored clips, 14 colored rings, 6 brown base rings and IFU)	1 lead set	TC70: 860315, 860352, 860353; TC50: 860310; TC30: 860306, 860354, 860355		•
989803151641	Complete lead set, IEC (2 arm leads, 2 leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 2 banana post adapters, 2 Trident yoke labels (C1-C3, C4-C6), 4 colored clips, 14 colored rings, 6 white base rings and IFU)	1 lead set		•	•
989803151651	Long complete lead set, AAMI (2 long arm leads, 2 long leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 2 Trident yoke labels (V1-V3, V4-V6), 4 colored clips, 14 colored rings, 6 brown base rings and IFU)	1 lead set			•
989803151661	Long complete lead set, IEC (2 long arm leads, 2 long leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 10 banana post adapters, 2 Trident yoke labels (C1-C3, C4-C6), 4 colored clips, 14 colored rings, 6 white base rings and IFU)	1 lead set		•	•
989803157481	DigiTrak XT patient cable – short	1 patient cable	DigiTrak XT: 860322	•	•
989803157491	DigiTrak XT patient cable – medium	1 patient cable		•	•
989803157501	DigiTrak XT patient cable – long	1 patient cable		•	•
989803175891	PageWriter TC20 10-lead patient cable IEC 3.4 m (11') long	1 patient cable	TC20: 860332 ST80i: 860343, 860351	•	•
989803175901	PageWriter TC20 10-lead patient cable AHA 3.4 m (11') long AAMI	1 patient cable			•
989803175911	PageWriter TC20 long 10-lead patient cable IEC 3.8 m (12.5') long	1 patient cable		•	•
989803175921	PageWriter TC20 long 10-lead patient cable AHA 3.8 m (12.5') long AAMI	1 patient cable			•
989803180121	ST80i Stress 12-lead lead set, grabber, AAMI	1 lead set			•

Ordering code	Description	Packaging	Use with	CE	FDA
989803180131	ST80i Stress 12-lead lead set, grabber, IEC	1 lead set	ST80i: 860343, 860351	•	
989803180141	ST80i Stress 12-lead lead set, snaps, AAMI	1 lead set			•
989803180151	ST80i Stress 12-lead lead set, snaps, IEC	1 lead set		•	
989803180161	ST80i Stress 12-lead lead set, grabber, AAMI long	1 lead set			•
989803180171	ST80i Stress 12-lead lead set, grabber, IEC long	1 lead set		•	
989803180181	ST80i Stress 12-lead lead set, snaps, AAMI long	1 lead set			•
989803180191	ST80i Stress 12-lead lead set, snaps, IEC long	1 lead set		•	
989803184921	Pagewriter TC10 standard 10-lead patient cable IEC 3.4 m (11') long	1 patient cable	TC10: 860392	•	
989803184931	Pagewriter TC10 standard 10-lead patient cable AHA 3.4 m (11') long AAMI	1 patient cable			
989803184941	Pagewriter TC10 long 10-lead patient cable IEC 3.8 m (12.5') long	1 patient cable		•	
989803184951	Pagewriter TC10 long 10-lead patient cable AHA 3.8 m (12.5') long AAMI	1 patient cable			

Patient cables



989803151631



989803175891



989803175901



989803175921_a



989803175921_b



989803180151_a



989803180151_b



989803180151_c



989803180191_a




989803180191_b



989803180191_c



989803184921



989803184941



Patient data cables

Ordering code	Description	Packaging	Use with	CE	FDA
989803158481	Class A USB patient data cable Length: 6.6' (2 m)	1 cable = 1 sales unit	TC70: 860315	•	•
989803164281	Class B USB patient data cable Length: 6.6' (2 m)	1 cable = 1 sales unit	TC70: 860315, 860316, 860352, 860353; TC50: 860310, 860311; TC30: 860306, 860307, 860354, 860355	•	•



989803164281

Replacement lead sets

Ordering code	Description	Packaging	Use with	CE	FDA
989803151671	Chest lead set, AAMI and IEC (1 x 1 into 3 Trident leads, 12 x 2 colored rings, 3 brown base rings, 3 white base rings (both AAMI and IEC), 3 banana post adapters, 4 Trident yoke lead description labels (V1-V3, V4-V6, C1-C3, C4-C6) and IFU)	1 lead set	TC70: 860315, 860352, 860353; TC50: 860310; TC30: 860306, 860354, 860355	•	•
989803151691	Long chest lead set, AAMI and IEC (1 x 1 into 3 Trident leads, 3 x 2 colored rings, 3 brown base rings, 3 white base rings (both AAMI and IEC), 3 banana post adapters, 1 lead separator, 4 Trident yoke lead description labels (V1-V3, V4-V6, C1-C3, C4-C6) and IFU)	1 lead set		•	•
989803151711	Limb lead set, AAMI and IEC (2 arm leads, 2 leg leads, 5 x 2 colored rings and 13 x 2 colored clips for AAMI and IEC, 4 banana post adapters and IFU)	1 lead set		•	•
989803151731	Long limb lead set, AAMI and IEC (2 long arm leads, 2 long leg leads, 13 x 2 colored rings and 5 x 2 colored clips for AAMI and IEC, 4 banana post adapters and IFU)	1 lead set		•	•



989803151671



989803151711

Extended leads

Ordering code	Description	Packaging	Use with	CE	FDA
989803151701	16-lead spare parts kit (9 x 2 colored marking clips, 4 banana post adapters, 6 shorting plugs, 2 lead separators 2.54 DIA, 2 lead separators 3.4 DIA and IFU)	1 kit	TC70: 860315, 860352, 860353; TC50: 860310		
989803151751	Upgrade Kit 12-15/16 leads, AAMI (4 additional lead assembly, 6 shorting plugs, 2 lead separators, 6 colored clips, 6 colored rings, 4 brown base rings and IFU)	1 set			•
989803151761	Upgrade Kit 12-15/16 leads, IEC (4 additional lead assembly, 6 shorting plugs, 2 lead separators, 4 banana post adapters, 6 colored clips, 6 colored rings, 4 white base rings and IFU)	1 set		•	•
989803151771	Upgrade Kit 12-15/16 long leads, AAMI (4 additional long lead assembly, 6 shorting plugs, lead separators, 6 colored clips, 6 colored rings, 4 brown base rings and IFU)	1 set			•
989803151781	Upgrade Kit 12-15/16 long leads, IEC (4 additional long lead assembly, 6 shorting plugs, 2 lead separators, 4 banana post adapters, 6 colored clips, 6 colored rings, 4 white base rings and IFU)	1 set		•	•



989803151701a

989803151701b

Electrode connectors

Ordering code	Description	Packaging	Use with	CE	FDA
989803129241	Alligator clips, IEC For tab electrodes	10 adapters	TC70: 860315, 860352, 860353; TC50: 860310; TC30: 860306, 860354, 860355;	•	•
989803166031	Clear tab/snap adapter Enables the connection from 3–4 mm post banana plug coming from the cardiograph to tab and snap style electrodes	10 adapters	TC20: 860332; TC10: 860392	•	•



989803166031a




989803166031b



989803166031c

Disposable ECG electrodes


Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803100441	13943B	Adult solid gel tab electrode for resting ECG 1-3/8 x 7/8" (35 x 22 mm)	100 per pack; 10 packs per case = 1,000 electrodes	TC70: 860315, 860352, 860353; TC50: 860310; TC30: 860306, 860354, 860355; TC20: 860332; TC10: 860392 Clinical research: 860297	•	•
989803101301	40420A	High-tack foam, wet gel ECG electrode Used for diagnostic, stress, resting, and monitoring ECG 1.1 x 2.2" (28 x 55 mm)	30 per pack; 10 packs per case; = 300 electrodes	ST80i: 860343, 860351 TC10: 860392	•	•
989803129091	M4612A	Solid gel snap ECG electrode	5 per pack; 60 packs per case; = 300 electrodes	DigiTrak XT: 860322	•	•
989803129101	M4613A	Solid gel snap electrode	30 electrodes per pouch; 10 pouches per case		•	•
989803192551		Clear tape wet gel ECG electrode for monitoring adult ECG 50 mm (2") round	30 per pack; 10 packs per case; = 300 electrodes	StressVue: 860295; DigiTrak XT: 860322	•	•



989803101301



989803129101a



989803134771



989803192541



989803192551

Limb plate electrodes

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE FDA
989803101691	40494E	Limb clamp electrode, clothespin type for 1/8" post end leads or 4 mm banana plug leads, screw connector	4 electrodes	All cardiographs with standard 1/8" post end or 4 mm banana plug leads	• •



989803101691

Recording paper

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803104771	M1708A	PageWriter XL thermal paper Z-fold (216 x 279 mm) with red grid and no header Size: Letter	1 case = 8 pads = 1,600 sheets	ST80i: 860343, 860351		•
989803104791	M1710A	PageWriter XL thermal paper Z-fold (210 x 300 mm) with red grid and no header Size: A4	1 case = 8 pads = 1,600 sheets		•	•
989803106261	M2481A	Thermal paper for PageWriter cardiographs Z-fold with 206 mm light red grid and 2 7/8" (7.3 cm) header Size: Letter	1 case = 8 pads = 1,600 sheets	TC70: 860315, 860352, 860353; TC50: 860310; TC30: 860306, 860354, 860355; ST80i: 860343, 860351	•	•
989803106271	M2483A	Thermal paper for PageWriter cardiographs Z-fold with 233 mm light red grid and 6.7 cm header Size: A4	1 case = 8 pads = 1,600 sheets		•	•
989803106281	M2485A	PTP brand anti-fade thermal paper for PageWriter cardiographs Z-fold with 206 mm light red grid and 2 7/8" header (7.3 cm) Size: Letter	1 case = 8 pads = 1,600 sheets		•	•
989803136421	M3707A	Thermal paper for PageWriter Trim Z-fold with light red grid and 2 7/8" (7.3 cm) header Size: Letter	1 case = 16 pads = 1,600 sheets		•	•
989803136431	M3708A	Thermal paper for PageWriter Trim Z-fold with light red grid and 7.3 cm header Size: A4	1 case = 16 pads = 1,600 sheets		•	•
989803185271		Thermal roll paper, 110 mm wide	1 case = 10 rolls	TC10: 860392	•	

Pouches and cases

Ordering code	Description	Packaging	Use with	CE	FDA
989803136651	DigiTrak Plus pouch	1 pouch	DigiTrak Plus: 860277	•	•
989803137971	StressVue velcro belt with pouch for wired patient module	1 pouch	StressVue wired patient module	•	•
989803144351	StressVue pouch for exchange wireless patient module labeled "X12+"	1 pouch	StressVue wireless patient module labeled X12+	•	•
989803153451	DigiTrak XT pouch	1 pouch	DigiTrak XT: 860322	•	•
989803158191	DigiTrak XT belt clip	1 each		•	•
989803158201	DigiTrak XT belt clip - 10 pack	10 per pack		•	•
989803166501	PageWriter TC30 and TC50 Cardiograph keyboard membrane cover	5 per pack	TC50: 860310; TC30: 860306, 860354, 860355; TC20: 860332		
989803166511	PageWriter TC70 Cardiograph keyboard membrane cover	5 per pack	TC70: 860315, 860352, 860353		
989803180801	Stress patient belt to secure patient interface module (PIM) for stress system (adjustable size)	1 each	ST80i: 860343, 860351	•	•
989803180811	Stress PIM holder	1 each	ST801: 860343, 860351, 453564277261, 453564277271	•	•



989803153451



989803166511

Cardiograph batteries

Ordering code	Description	Packaging	Use with	CE	FDA
989803170371	Single battery PW-TC series Minimum shelf life: 6 months	1 battery	TC70: 860315, 860352, 860353; TC50: 860310; TC30: 860306, 860354, 860355		
989803185291	Pagewriter TC10 battery	1 battery	TC10: 860392		
989803193431	Lithium ion battery 7.252 Wh	1 battery	ECG100		



989803185291

Imaging supplies

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803101601	40483A	Aquasonic gel (not available in China)	12 bottles (8.5 oz/250 ml)	All ultrasound imaging, Avalon US Transducer: M2736A; Avalon Wireless Transducer: M2720A	•	•



989803101601

TEE supplies


Ordering code	Alternative ordering code	Description	Packaging	Use with	CE FDA
989803105311	M1828A	Disposable TEE bite guard, adult	24 guards	TEE Transducers: 21362A, 21362B, 21367A, 21369A; HP ultrasound probe: 21362C; Biplane transducer: 21363A; TEE probe: 21364A; MPT transducer: 21478A; Sonos 2000C: M2406A; Omniplane III TEE probes: 21378A, 21478A	
989803133811	M4604A	Tip protector for Omniplane III	24 tip protectors	Omniplane III TEE probes: 21378A, 21478A	•



989803101621



989803106111



989803133811

List of host products

Ordering code	Alternative ordering code	Description
860277		DigiTrak Plus holter recorder
860295		StressVue exercise stress testing system
860297		Clinical research cardiograph
860306		PageWriter TC30 cardiograph
860307		PageWriter TC30 upgrade
860310		PageWriter TC50 cardiograph
860311		PageWriter TC50 upgrade
860315		PageWriter TC70 cardiograph
860316		PageWriter TC70 upgrade
860322		DigiTrak XT holter recorder
860332		PageWriter TC20 cardiograph
860333		PageWriter TC20 upgrade
860343		ST80i Stress Test System
860351		ST80i Stress Test System upgrade
860352		PageWriter TC70 government bundle with trolley
860353		PageWriter TC70 government bundle without trolley
860354		TC30 government bundle with trolley
860355		TC30 government bundle without trolley
860392		PageWriter TC10 cardiograph
860393		PageWriter TC10 upgrade

Ordering code	Alternative ordering code	Description
	21362A	TEE transducer
	21362B	TEE transducer
	21362C	HP ultrasound transducer
	21363A	Biplane transducer
	21364A	TEE transducer
	21367A	TEE transducer
	21369A	TEE transducer
	21378A	Omniplane III TEE transducer
	21478A	Omniplane III TEE transducer
	M2406A	Sonos 2000C ultrasound
862491	M2720A	Avalon CTS cordless fetal transducer system
989803143691	M2736A	Avalon US transducer

