







www.boehm-elektromedizin.de Tel: 030 408 999 0







15:2 mode 30:2 mode Compression Mode



## Battery

E6 features two durable batteries with 90mins working time.



## EtCO<sub>2</sub>

EtCO<sub>2</sub> module for monitoring the patients' carbon dioxide concentration.



#### Portable Carrying Bag

Machine net weight 8.9kg, Compact and lightweight.





#### Easy to install

- 1. EtCO2 port, charging port, USB port. 2. Battery charging dock.
- 3. Patient fixing strap.
- 4. Two compression plates included, 15mm and 25mm for different compression depth. 5. Portable carrying bag, easy to pack.


























### Easy to operate

- 1. Put the back plate underneath the patient.
- 2. Attach the support legs to the back plate and lock against each other.
- 3. Automatically positioning on the chest, one-key compression.
- 4. Tighten the fixing strap.








#### Touchscreen

- 3.5" color touchscreen.
- Simple Operation, Clarity and Intuitive User Navigation.



Ambulanc is one of the leading global manufacturer of Emergency and Transport ventilators. We are dedicated to innovation in the fields of emergency and hospital transportation for ventilation.

★ Headquarter

Headquartered in Shenzhen, China, Ambulanc possesses a sound distribution and service network with our distributors in more than 60 countries in North and Latin America, Europe, Africa and Asia-Pacific. While improving the quality of the products and service, we help in reducing its cost to save human lifes in emergency.

Since its foundation in 2001, Ambulanc's development has been driven by innovation with strong R&D inspired by the needs of our customers, we adopt advanced technologies and transform them into accessible innovation, bringing emergency solution for peoples.

# Ambulanc(shenzhen)Tech.Co.,Ltd. € 🚾



Sales Countries

Add: 3th Floor, Block C, Building #5, Skyworth Innovation Industry Park, Tang Tou 1st Road, Shiyan Town, Baoan District, Shenzhen 518108, China.

Tel: +86-755 26073861 Fax: +86-755 23016012

Web site: www.ambulgroup.com E-mail:info@ambu-lanc.com